

Caltest Container Request Form

Bottle Order # _____

FORM FILLED OUT BY: _____

DATE REQUESTED: _____

REQUIRED RECEIPT DATE: _____

Client: _____
 Address: _____
 Delivery Attn: _____

CLIENT PICK UP TIME _____
 CALTEST TO DELIVER _____
 SHIP Ground OverNight

QTY

PLASTIC CONTAINERS

- Metals, Pint (~500mL) with **HNO₃** preservative
- Metals, Dissolved, Pint (~500mL) unpreserved
- Cyanide, TOT or WAD, **Amber** Pint (~500mL) with **NaOH** pres. Sodium Thiosulfate dechlorination kit
- Sulfamic acid kit for NO₃/NO₂ pretreatment (incl. 3x NaOH)
- Pint (~500 mL) with **H₂SO₄** preservative for NH₃, TKN, Total PO₄, NO₃-NO₂-N, COD
- Half Pint (~250 mL) with **H₂SO₄** preservative for one or two of the above
- Coliform (~150 mL) sterile with **Na₂S₂O₃** dechlorination for P/A, MPN, E.COLI or SPC/HPC Bacteria Analysis
- Coliform (~290 mL) sterile with **Na₂S₂O₃** dechlorination for MPN or E.COLI requiring analysis in duplicate
- 549 Diquat (Title 22) – **Plastic Amber** Liter with **Na₂S₂O₃** for dechlorination
- Hexavalent Chromium (CrVI) by EPA 218.6 (for Drinking Water or Region 5) - 60mL plastic w/ **NH₄SO₄/NH₄OH** buffer
- Chlorophyll A (Algal Biomass) **Plastic Amber** Liter unpreserved
- Cubitainer for % Survival Bioassays: 5 gallon 2.5 gallon 1 gallon
- Gross Alpha Gross Beta Ra226 Ra228 Uranium
- Bromate - Plastic 50mL with **EDA Solution** preservative
- Unpreserved Plastic 1/2Pint(~250mL) Quart (~950mL) Half Gallon
- Pint (~500mL) Liter (1000mL)* *required for Lead&Copper Rule and TSS/TDS/TS

FILLED BY: _____

CHECKED BY: _____

DELIVERED BY: _____

GLASS CONTAINERS

- Phenol by EPA 420.1 - 250 mL amber glass with **H₂SO₄** preservative
- 1631 Mercury Ultra Trace Level – 500mL glass double bagged Tot (HCl pres) Dis (unpres)
- 1630 Methyl Mercury – 500mL **Amber** glass double bagged Tot (HCl pres) Dis (unpres)
- 1630 Travel Blank for Methyl Mercury – 500mL **Amber** glass with **HCl** preservative and double bagged (keep cool)
- Mercury free water for 1630 or 1631 Field Blank – 500mL glass UNPRESERVED and double bagged
- Amber Liter (AL) wide mouth for Semi-Volatile Organics – for the following method * **Na₂S₂O₃** for dechlorination if needed
- 3 AL per analysis:** 625.1-BNA's*, PAHs*, 625.1-OCP's* (formerly 608), 625.1-OPP's* (formerly 614), 625.1/8270 Pyrethroids*, 632.
- Reduced AL Qty may be acceptable depending on sample/analysis count; please inquire. **2 AL per analysis:** Diesel/MotorOil
- Oil & Grease by EPA 1664, Amber Liter (1000mL) with **HCl** preservative
- Oil & Grease by EPA 1664, with **HCl**, Industrial Waste matrix or Compositing 500mL 250mL 125mL 60mL
- General Physical (Title 22) – Amber Liter (1000mL) unpreserved
- 1613 Dioxins/Furans - Amber Liter (Collect in Duplicate) **Na₂S₂O₃** for dechlorination if needed
- Tributyltin Amber Liter 1000mL with **HCl** preservative (Collect in Duplicate)
- 525 or 507 (525 A&S only) Amber Liter 1000mL with **Na₂S₂O₃** for dechlorination and vial of **HCl** for preservative (Collect in Duplicate)
- 547 - 125mL Amber Glass or 40mL Amber VOAs; both with **Na₂S₂O₃**
- 548 - Amber 250mL with **Na₂S₂O₃** for dechlorination (Collect in Duplicate)
- Sulfide - 250mL clear glass bottle preserved for
- Soil Jars: 4 oz. 8 oz. 32 oz. Total Poly Jar Dissolved Amber 8 oz.
- Sulfite (SO₃) – 500 mL Amber Glass with **EDTA** solution for preservative
- Dissolved Organic Carbon (DOC) - Amber Glass 125mL unpreserved

GLASS SEPTA BOTTLES, 40 mL VIALS

- 515 - (4) 60mL Amber VOAs with **NaSulfite** for dechlorination or Amber 125mL Glass with **Sulfite** for dechlorination
- Triplicate 40mL VOAs with **HCl** preservative for: 8260, MTBE, BTEX, Gas, 524.2 524.2 **Ascorbic Acid** dechlorination kit
- 524.2 THMs Only - Triplicate 40mL VOAs with **Na₂S₂O₃** for dechlorination
- 624.1- Two VOAs preserved with HCl, Two VOAs unpreserved, **total of 4 VOAs.**
- 624.1 - Two VOAs preserved with HCl, Two VOAs unpreserved, and Two VOAs client preserved to pH 4-5 for a **total of 6 VOAs.**
- *624 Dechlorination kit (250 mL glass) with **Sodium Thiosulfate** if sample is chlorinated (**Na₂S₂O₃** for 624)
- Travel Blank, triplicate preserved 40mL clear VOA with **HCl** preservative.
- EDB/DBCP: (3) Amber VOAs with **NH₄Cl** preservative + **Phosphate Buffer** (504/551) or Amber VOA with **Ascorbic/Maleic acid** (524.3)
- 505 or 508 - (4) Amber 40mL vials (505) or Duplicate Amber Glass Liters (508); both with **Na₂S₂O₃** for dechlorination
- 531 - Two Amber VOAs preserved with **KH₂ Citrate** and **Na₂S₂O₃** for dechlorination
- 552 Haloacetic Acids (HAA) – Triplicate Amber Glass VOA vials or 125mL amber glass with **NH₄Cl** preservative
- Formaldehyde - Triplicate Amber 40mL VOAs preserved with **NH₄Cl/CuSO₄**
- Total Organic Carbon (TOC) – Triplicate Amber Glass 40mL VOA set with **HCl** preservative

OTHER ITEMS

- Filter Kits for: Nutrients/Merc ICPMS Metals
- Cubitainers filled with DI Water 5 gallon 2.5 gallon 1 gallon
- Ice Chest(s) Container Labels Chain of Custody Forms Customized labels per attached COC
- Other: _____